


Various Methods to Access PACS

With CST Cerner live at BC Children's and Women's Hospitals, there are now several different ways to access PACS to view Medical Imaging exams

1. Via the PHSA Philips PACS icon on the desktop

- Launch PHSA Philips PACS using the desktop icon



- CST Cerner patients are listed under the Cerner-PHSA organization in Philips iSite
- This method only works if you have access to Philips iSite

2. Via Powerchart

- Load the patient's chart in Powerchart
- Launch Philips iSite PACS using the PHSA PACS button in the toolbar



- Philips iSite PACS will launch in context with the patient
- This method only works if the patient already had a medical imaging order in CST Cerner
- If the patient does not have any medical imaging orders in CST Cerner yet, then the user needs to use the legacy Cerner Millennium to launch iSite

3. Via Results Review

- In Powerchart, navigate to results review
- Open the exam in question
- Launch Philips iSite PACS using the "view image" button in the toolbar



- Philips iSite PACS will launch in context with the patient and the exam in question
- This method only works if the patient already had a medical imaging order in CST Cerner and is meant for opening the images associated to that order
- If the patient does not have any medical imaging orders in CST Cerner yet, then the user needs to use the legacy Cerner Millennium to launch iSite

4. Via legacy Cerner Millennium

- Open legacy Cerner Millennium Powerchart



- Load the patient's chart
- Launch Philips iSite PACS using the iSite button in the toolbar

