

HEMOCHRON SIGNATURE ELITE (ACT+) QUALITY CONTROL RECORD FORM

Notes:

- 1. Start new sheet with new lot of cuvettes or new lot of QC materials.
- 2. Record lot numbers and expiry dates of cuvettes and QC materials.
- 3. Perform 2 levels of liquid QC every 30 Days.
- 4. Perform 2 levels of liquid QC with new shipment and/or new lot # of cuvettes

Hemochron Jr ACT+ Lot# Fridge Exp: Delivery Date QC PASS Y/N		INSTRUMENT CHECKS						LEVEL 1 QC (Normal)				LEVEL 2 QC (Abnormal)				OTHER	
		Instrument 1: SN11734 Instrument 2: SN11749						Lot#				Lot#					
								Exp:				Exp					
								Range				Range					
Date	Time	Instrument 1			Instrument 2			Instrument	QC PASS	Instrument	QC PASS	Instrument	QC PASS	Instrument	QC PASS	Patient Test	Operator ID
		Battery	Temp	EQC	Battery	Temp	EQC										