

## Urine Dipstick Emergency Dept: Form D - Patient Result Record Form

### Urinalysis C&W Chemstrip 10A

**Patient Name:** \_\_\_\_\_ **Date:** \_\_\_\_\_  
**Time Collect:** \_\_\_\_\_ **Time Test:** \_\_\_\_\_  
**Volume:** \_\_\_\_\_ mLs

**Circle: specimen type/color/clarity**

**Specimen Type:** midstream- bagged – catheter - other  
**Color:** yellow – orange – red – green – blue – brown - other  
**Clarity:** clear – cloudy - slightly cloudy – turbid - other

Read @ (seconds)	Compare reagent area on test strip at specified times to color chart on test strip bottle. Circle result.							
<b>SG</b>	60	1.000	1.005	1.010	1.015	1.020	1.025	1.030
<b>pH</b>	60			5	6	7	8	9
<b>Leukocytes</b>	60-120	Leucocytes/uL		neg	1+ 10-25	2+ 75	3+ 500	
<b>Nitrite</b>	60			neg	pos			
<b>Protein</b>	60	g/L		neg	1+ 0.3	2+ 1.0	3+ 5.0	
<b>Glucose</b>	60	mmol/L		normal	1+ 2.8	2+ 5.5	3+ 17	4+ 56
<b>Ketones</b>	60	mmol/L		neg	1+ 1.5	2+ 5	3+ 15	
<b>Urobilinogen</b>	60	umol/L		normal	1+ 17	2+ 70	3+ 140	4+ 200
<b>Bilirubin</b>	60	umol/L		neg	1+ 17	2+ 50	3+ 100	
<b>Blood</b>	60	Erythrocytes/uL		neg	1+ 5-10	2+ 25	3+ 50	4+ 250
				intact	1+ 10	2+ 25	3+ 50	4+ 250
				lyzed	1+ 10	2+ 25	3+ 50	4+ 250

### Urinalysis C&W Chemstrip 10A

**Patient Name:** \_\_\_\_\_ **Date:** \_\_\_\_\_  
**Time Collect:** \_\_\_\_\_ **Time Test:** \_\_\_\_\_  
**Volume:** \_\_\_\_\_ mLs

**Circle: specimen type/color/clarity**

**Specimen Type:** midstream- bagged – catheter - other  
**Color:** yellow – orange – red – green – blue – brown - other  
**Clarity:** clear – cloudy - slightly cloudy – turbid - other

Read @ (seconds)	Compare reagent area on test strip at specified times to color chart on test strip bottle. Circle result.							
<b>SG</b>	60	1.000	1.005	1.010	1.015	1.020	1.025	1.030
<b>pH</b>	60			5	6	7	8	9
<b>Leukocytes</b>	60-120	Leucocytes/uL		neg	1+ 10-25	2+ 75	3+ 500	
<b>Nitrite</b>	60			neg	pos			
<b>Protein</b>	60	g/L		neg	1+ 0.3	2+ 1.0	3+ 5.0	
<b>Glucose</b>	60	mmol/L		normal	1+ 2.8	2+ 5.5	3+ 17	4+ 56
<b>Ketones</b>	60	mmol/L		neg	1+ 1.5	2+ 5	3+ 15	
<b>Urobilinogen</b>	60	umol/L		normal	1+ 17	2+ 70	3+ 140	4+ 200
<b>Bilirubin</b>	60	umol/L		neg	1+ 17	2+ 50	3+ 100	
<b>Blood</b>	60	Erythrocytes/uL		neg	1+ 5-10	2+ 25	3+ 50	4+ 250
				intact	1+ 10	2+ 25	3+ 50	4+ 250
				lyzed	1+ 10	2+ 25	3+ 50	4+ 250

## REVISION & APPROVAL LOG

Version	Revision Type	Description of Change	Revision Date	Technical Approval	Medical Approval
1.0		New document	Dec 2013	Elvira Kozak	Dr. Cathy Halstead
1.1	Minor	Revisions	19 May 2016	Elvira Kozak	
1.1	Minor	Document title and number change. Upload to QMS document control	28 Dec 2016		Dr. Benjamin Jung

**Attention: This document is published on the BCCW SharePoint website**

**Revisions made to this document require an update to the corresponding document published on BCCW SharePoint website**