





Site Applicability


BC Women's Hospital – Labour and Delivery, Delivery Rooms

Set Up

Remove:

- Ear probe and tips
- EZ scrub packages
- Extra thermometer covers




Area	Instructions	Pictures
Sink area	1 box gloves for primary RN size	
FHOL warmer	Leave in hall until needed for delivery (NICU to decide supplies in drawers)	
"Covid+ /PUI Delivery Cart"	Leave in hall until needed for delivery.	<p>Delivery cart: Reduced items for Covid + /PUI</p>
"Covid +/PUI Nursing Cart"	In room. Drawer 3 C/S supplies placed on top of cart in hall and bring if needed.	<p>Nursing cart: Reduced items for Covid + /PUI</p>
Anesthesia cart	Leave in hall or ante-room for possible epidural use- anesthesia to collect equipment from cart outside room for pt epidural	
Ante-room	Keep PPE stock as Labeled - No extra supplies	


Baby cot	Reduce supplies as needed	
Head of Bed	Keep 1 of each emergency equipment pieces as labeled	

Revised Isolation Delivery Cart

**Delivery cart:
Reduced items
for
Covid + /PUI**

Items removed- pudendal block, blood collection for studies, vacuum and forceps and under buttocks drape, extra drape. Extra items to be stored in regular delivery cart outside room




Drawer	Contents	Picture
1	1 pack sponge, 1 plastic placenta bag, 1 pair goggles, 1 instrument count sheet	
2	2 blood gases and needles with bags, 1 yellow clamp, 2 sets of 2-0 and 3-0 sutures, 2 1% lidocaine vials and 1 syringe and 22 gauge needle	
3	2 size each glove- 5.5, 6, 6.5, 7, 7.5, 8	




4	Delivery set	
5	OB pack and single gown	

Revised Isolation Nursing Cart

**Nursing cart:
 Reduced items
 for
 Covid + /PUI**

Items removed- IUPC

Drawer	Contents	Pictures
Top	C/S drawer out of room bring in if needed for prep (in plastic bag)	
1	4 save a day tray 3 wash cloths, 1 pair goggles, 1 towel, 1 gown, 1 box Kleenex, 1 package wipes	
2	2 each of sterile gloves S,M,L , 2- 3ml 22 gauge syringe with needle, 1 spare box thermometer probe 1 ml syringe with needle, 2 blunt filter needles, 2- 5ml syringe, 1 ampule breaker, 4 med labels, 6 lubricant, 1-oxytocin 50ml, 1 measuring tape, 4- oxytocin 10ml, 1 vit K, 10 alcohol swabs 1 erythro, 2 blunt needles, 1 amnio hook,	

<p>3</p>	<p>1 pack blue pads, 1 pack maternity pads, 1 disposable undies, 1 cavi-wipes</p>	
<p>4</p>	<p>1 soaker pad, 1 catheter tray, 1 foley catheter, 1 I&O, 1 leg cath sticker</p>	
<p>5</p>	<p>1 FECG electrode attachment, 1 monitoring belt, 1 monitoring cables for ultrasound, toco, IUPC, FECG MicroClave extension set, IV supplies- 1L N/S, roll tape, mainline tubing, 2 Alcohol swab, #18 GA Insyte, Gauze, bloodwork tubes- 2 purple, green, blue, gold, tourniquet, vacutainer adaptor, black and blue</p>	

Ante Room Set Up

Ante room set up in L&D and UCC are standardized to be the same.

All ante rooms have shelves that are labelled for supplies.

No other equipment or supplies are to be in the ante rooms other than what is labelled.



Version History

DATE	DOCUMENT NUMBER and TITLE	ACTION TAKEN
09-Apr-2020	C-06-15-60607 Delivery Room Set Up: COVID-19 Positive Or PUI Patient	Developed by BCW Maternal Newborn Program; Approved by Professional Practice Director

Disclaimer

This document is intended for use within BC Children's and BC Women's Hospitals only. Any other use or reliance is at your sole risk. The content does not constitute and is not in substitution of professional medical advice. Provincial Health Services Authority (PHSA) assumes no liability arising from use or reliance on this document. This document is protected by copyright and may only be reprinted in whole or in part with the prior written approval of PHSA.